



## **XML ECG Anonymous File Converter** (available in English only)

**Click on a blue button to select a file**

---

**This software is provided as-is, without any express or implied warranty. In no event shall Philips Medical Systems be held liable for any damages arising from the use of this software.**

---

**The XML ECG Anonymous File Converter application removes all personal patient identification information and institution and other location information from an ECG file. The resulting generic ECG file can then be used for research or other purposes. The application also converts any standard Philips ECG XML file with UTF-16 character encoding to an XML ECG file with UTF-8 character encoding. This conversion process reduces the overall size of the XML ECG file by half.**

**This converter application has been designed for and exclusively tested with US English ECG files. This application should not be used with ECG files that require full UTF-16 encoding.**

**This application is compatible with XML ECGs generated by the Philips PageWriter Trim and PageWriterTouch cardiographs, and the Philips Intellivue Bedside Monitor and Information System (with 12-lead capability), or with XML ECGs exported by the Philips TraceMasterVue ECG Management System.**

**The XML ECG Anonymous File Converter application operates with Microsoft XML Core Services (MSXML 4.0 library) version 4.0 or higher. Click the button below for more information.**

**More  
Info**

**Click to read the program documentation as a Readme.txt file**

**More  
Info**

**Click to run the setup.exe installer**

**More  
Info**

**Click to read the program documentation for the Microsoft  
MSXML 4.0 Service Pack 2 program**

---